


96-well Deep Well Microplate

Critical Dimensions

Drg. 201116-B / Issue 2



<i>L</i>	128.15	± 0.10
<i>W</i>	85.85	± 0.10
<i>H</i>	44.40	± 0.10
<i>Hs</i>	4.00	± 0.10
<i>X1</i>	14.00	± 0.10
<i>Xn</i>	113.50	± 0.10
<i>Y1</i>	11.00	± 0.10
<i>Yn</i>	74.50	± 0.10

<i>P</i>	9.00	± 0.05
<i>D</i>	8.25	± 0.05
<i>H350</i>	15.30	± 0.10
<i>H1</i>	24.50	± 0.10
<i>H2</i>	39.90	± 0.10
<i>C</i>	12	
<i>R</i>	8	

<i>Volume 350µl</i>	750µl
<i>Volume 1ml</i>	1400µl
<i>Volume 2ml</i>	2500µl
<i>Well form</i>	V
<i>Material</i>	PP
<i>Cv Max</i>	0.50
<i>Cw Max, Cl Max</i>	1.00

Dimensions in mm