


Aqueous SEC 2 Std. (AL0-3043)



Column: PolySep-SEC 0 μ m Polymeric (synthetic) - Acrylic (unmodified) SEC 1000

Dimensions: 300 x 7.8 mm ID

Order No:

Elution Type: Isocratic

Eluent A: water

Pct. A: 100

Flow Rate: 0.8 mL/min

Temperature: 22 °C

Detection: Refractive Index @ nm

SAMPLE ANALYTE

1 ethylene glycol

SPE:

For more information call your Phenomenex Representative at (310) 212-0555

Copyright 2000 Phenomenex Inc.

Ref. No.: 14753/QUA/173129-3