


Technique:  
TN-2014

Application Note:

## Non-Halogenated Organics using GC-FID


### ZB-WAX 30m x 0.53mm x 1.0µm



### Compound List

- |                            |                                 |
|----------------------------|---------------------------------|
| 1. Methanol (Solvent)      | 8. Isobutanol                   |
| 2. Ethanol                 | 9. 1,4-Dioxane                  |
| 3. Acetonitrile            | 10. Methylmethacrylate          |
| 4. Diethyl ether           | 11. 4-Methyl-2-pentanone (MIBK) |
| 5. Propanitrile            | 12. 2,2-Dimethoxybutane         |
| 6. 2-Butanone (MEK)        | 13. Ethylmethacrylate           |
| 7. 2-Methyl-2-propanitrile | 14. Acrylamide                  |

### ZB-5 30m x 0.53mm x 1.5µm



www.phenomenex.com  
 Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department by telephone, fax or e-mail: international@phenomenex.com.

**phenomenex®**  
 ...breaking with tradition

|                                                                                         |                                                                               |                                                                          |                                                                                |                                                                             |                                                                         |                                                                            |
|-----------------------------------------------------------------------------------------|-------------------------------------------------------------------------------|--------------------------------------------------------------------------|--------------------------------------------------------------------------------|-----------------------------------------------------------------------------|-------------------------------------------------------------------------|----------------------------------------------------------------------------|
| <b>USA</b><br>tel.: (310) 212-0555<br>fax: (310) 328-7768<br>email: info@phenomenex.com | <b>Puerto Rico</b><br>(800) 541-HPLC<br>(310) 328-7768<br>info@phenomenex.com | <b>Canada</b><br>(800) 543-3681<br>(310) 328-7768<br>info@phenomenex.com | <b>United Kingdom</b><br>01625-501367<br>01625-501796<br>ukinfo@phenomenex.com | <b>Germany</b><br>06021-58830-0<br>06021-58830-11<br>anfrage@phenomenex.com | <b>New Zealand</b><br>09-4780951<br>09-4780952<br>info@phenomenex.co.nz | <b>Australia</b><br>1800-553-929<br>1800-553-923<br>info@phenomenex.com.au |
|-----------------------------------------------------------------------------------------|-------------------------------------------------------------------------------|--------------------------------------------------------------------------|--------------------------------------------------------------------------------|-----------------------------------------------------------------------------|-------------------------------------------------------------------------|----------------------------------------------------------------------------|